



Programming: SimDesigner

Upon full participation and completion of this course, participants will be able to:

- Start up SimDesigner.
- Create learner events in a resource list.
- Create a scenario containing multiple frames.
- Create and utilize either a trend or transition time within a scenario.
- Convert existing legacy scenarios using SimDesigner.

Course Agenda

- | | |
|--|---|
| ① Welcome and Introduction | ⑧ Questions & Closing Remarks |
| ② Introduction to Scenario Development | ⑨ Adding Timers |
| ③ Start Up Procedure and User Interface Overview | ⑩ Automatic Responses |
| ④ Scenario Information | ⑪ Multiple Phase Considerations |
| ⑤ Phases and Resources | ⑫ Special Circumstances for Events |
| ⑥ Patient States | ⑬ User Defined Variables |
| ⑦ Scenario Conversion | ⑭ Implementing a Scoring and Grading System |

Educational Representative Responsibility

Your Laerdal® Educational Representative has multiple responsibilities during the introductory course:

- Help you understand Laerdal resources.
- Conduct an interactive agenda to optimize your understanding of the simulator system and facilitate the development of your proficiency in its operation.
- Facilitate discussions about how the simulator system can meet your specific program goals.

